PPE - FACE MASK - KN95 & 3 PLY





TAYASH TRADE AND IMPEX PVT LTD

Contents

- PPE
- KN95 Mask
- 3 PLY Mask



PPE – PERSONAL PROTECTION EQUIPMENT



Name of product	Materials	Color	Weight	Characteristics
Disposable medical isolation clothing	Non-woven fabric (SMS)	White	60g/m2	Soft touch, anti-static, waterproof, breathable, no metal, no latex, no penetration, effective isolation of microorganisms

Disposable medical isolation clothing conforms to GB19082–2009, chest opening, integral structure, can cover all clothing and exposed skin. Isolation clothing has the following characteristics:


1. is widely used in exposure to patients with contact-borne infectious diseases and their surroundings, such as VRE (Enterococcus vancomycin-resistant), C. difficile (Clostridium difficile), NV (e.g. viruses) and other intestinal pathogens, RSV (respiratory syncytial virus), whether it is a suspected or confirmed infection or a rating.

TAYASH HEALTH CARE









Packing information:

60pcs/carton, carton size:490*590*360mm 16cartons/pallet(960pcs), entire size: 1.2*1.0*1.58mm

MOQ – 20,000 PIECES





KN95 Disposable face mask (Without Breather Valve)

KN95 Mask made of non-woven fabric, comes with activated carbon layer, which is non-irritating and non-allergenic to the skin.


M-shaped nose clip, dust mask is cut in 3D fits almost any size/shape face, 360° protection. 3D three-dimensional ventilation space, more comfortable breathing.

Adopt electrostatic super filter element can filter out dust, germs, and particles in the air, Prevent the spread of influenza, respiratory system and other infectious diseases.

The filtration efficiency (dust-removing rate) of KN95 dust mask is more than 95% for non-oily particulate matter above 0.075 microns.

This respirator is made for 95% filtration efficiency against certain non-oil-based particles and aerosols.





- Execution Standard: GB2626-2006 KN95
- Filtration Efficiency: ≥95%.

Important:

Do not use when in places where the oxygen density is below 18%. -

Do not use if the mask is dirty and polluted

Do not wash and reuse

Do not use for more than 1 day

Storage

Storage Conditions: Storage temperature; -20 degrees (-4 degrees F) to

+40 degrees (+104 degrees F)

Storage humidity: less than 80% RH.

KN95 MASK - PACKING INFORMATION

	Packaging:	20 pcs /color box, 900 pcs / carton
	Box size:	14.8 * 7 * 13.8 cm
	MEAS SIZE:	60*40*50
Package	G.W./N.W. (carton)	9KGS/7.5 KGS
information	20ft container	201600 PCS
	40ft container	414000 PCS
	40hq HQ container	468000 PCS
	MOQ	100000 PCS





Medical Face Masks Ear Loop 3 ply

Face Mask, 3 ply is made of Polypropylene Spunbonded and Melt blown nonwoven fabrics with adjustable nose strip and flat or round elastic ear - loop to give basic protection to the user.

It's mainly used in cleanroom, critical environment, food handling/production and industrial manufacturing.



Components:


1st layer : Polypropylene Spunbonded nonwoven

2nd layer : Melt blown Nonwoven

3rd layer : Polypropylene Spun bond nonwoven

Nose strip: PE wire (Metal Free)

Ear-Loop: Flat elastic or round elastic (latex free)



Specifications:

Width across face: 175 mm + / - 2 mm

Standards:

Excellent breath ability - Non irritating, fibre glass free fabric - Non-woven spunbonded - 100% Polypropylene - Hypoallergenic - Splash resistant - Ultrasonic welding



Package information	Packaging:	20 pcs /color box, 900 pcs / carton	
	Box size:	14.8 * 7 * 13.8 cm	
	MEAS SIZE:	60*40*50	
	G.W./N.W. (carton)	9KGS/7.5 KGS	
	20ft container	201600 PCS	
	40ft container	414000 PCS	
	40hq HQ container	468000 PCS	
	MOQ	100000 PCS	

REMYA.S PH: +91-9946487980

tayashpurchase@gmail.com

TAYASH TAYASH TRADE AND IMPEX PVT LTD

Heathcare Add:No.51/649,K.T.Leyo Building,Vytilla-682019